

**Non-Human Primate
Pituitary Magnetic Bead
Panel 1**

96-Well Plate Assay

Cat. # NHPPT1MG-46K

MILLIPLEX[®] MAP

Non-Human Primate Pituitary Magnetic Bead Panel 1 96-Well Plate Assay

NHPPT1MG-46K

<u>TABLE OF CONTENTS</u>	<u>PAGE</u>
Introduction	2
Principle	3
Storage Conditions Upon Receipt	3
Reagents Supplied	4
Materials Required But Not Provided	5
Safety Precautions	5
Technical Guidelines	6
Sample Collection And Storage	7
Preparation of Reagents for Immunoassay	8
Immunoassay Procedure	10
Plate Washing	12
Equipment Settings	13
Quality Controls	15
Assay Characteristics	15
Troubleshooting Guide	17
Replacement Reagents	20
Ordering Information	21
Well Map	22

For Research Use Only. Not for Use in Diagnostic Procedures.

By purchasing this product, which contains fluorescently labeled microsphere beads authorized by Luminex[®] Corporation (“Luminex[®]”), you, the customer, acquire the right under Luminex[®] Corporation’s patent rights, if any, to use this product or any portion of this product, including without limitation the microsphere beads contained herein, only with Luminex[®] Corporation’s laser based fluorescent analytical test instrumentation marketed under the name of Luminex[®] 100[™] IS, Luminex[®] 200[™], Luminex[®] HTS, FLEXMAP 3D[®], MAGPIX[®].

INTRODUCTION

Pituitary hormones and other brain-derived proteins, such as hypothalamus neuropeptides, play very important roles in the regulation of various functions including metabolism, growth, and reproduction. Accurate measurement of these proteins to understand their new biological functions and molecular mechanisms of the functions are crucial. Traditional laboratory methods, such as RIA and ELISA are not able to measure multiple proteins with a small sample volume.

EMD Millipore recognizes the integral role that these proteins play and the MILLIPLEX[®] MAP Non-Human Primate Pituitary Magnetic Bead Panel 1 can be used for the simultaneous measurement of the following analytes in any combination: ACTH, AGRP, CNTF, FSH, GH, LH, and TSH. This multiplex assay can analyze these 7 proteins simultaneously and uses a small sample volume 25 μ L.

The MILLIPLEX[®] MAP Non-Human Primate Pituitary Magnetic Bead Panel 1 provides biomedical researchers with quality tools for the study of reproduction, growth, metabolic homeostasis, and pituitary-related diseases such as acromegaly, growth hormone deficiency, diabetes insipidus and pituitary tumors. MILLIPLEX[®] MAP enables you to investigate the modulation and expression of multiple analytes simultaneously, giving you the advantage of speed and sensitivity, and dramatically improving productivity.

EMD Millipore's MILLIPLEX[®] MAP Non-Human Primate Pituitary Magnetic Bead Panel 1 is the most versatile system available for metabolic hormone research.

- MILLIPLEX[®] MAP offers you the ability to:
 - Choose any combination of analytes from our panel of 7 analytes to design a custom kit that better meets your needs.
- A convenient "all-in-one" box format gives you the assurance that you will have all the necessary reagents you need to run your assay.

For Research Use Only. Not for Use in Diagnostic Procedures.

Please read entire protocol before use.

It is important to use same assay incubation conditions throughout your study.

PRINCIPLE

MILLIPLEX[®] MAP is based on the Luminex[®] xMAP[®] technology — one of the fastest growing and most respected multiplex technologies offering applications throughout the life-sciences, and is capable of performing a variety of bioassays including immunoassays on the surface of fluorescent-coded magnetic beads known as MagPlex[™]-C microspheres.

- Luminex[®] uses proprietary techniques to internally color-code microspheres with two fluorescent dyes. Through precise concentrations of these dyes, 100 distinctly colored bead sets can be created, each of which is coated with a specific capture antibody.
- After an analyte from a test sample is captured by the bead, a biotinylated detection antibody is introduced.
- The reaction mixture is then incubated with Streptavidin PE conjugate, the reporter molecule, to complete the reaction on the surface of each microsphere.
- The microspheres are allowed to pass rapidly through a laser, which excites the internal dyes marking the microsphere set. A second laser excites PE, the fluorescent dye on the reporter molecule.
- Finally, high-speed digital-signal processors identify each individual microsphere and quantify the result of its bioassay, based on fluorescent reporter signals.

The capability of adding multiple conjugated beads to each sample results in the ability to obtain multiple results from each sample. Open-architecture xMAP technology enables multiplexing of many types of bioassays reducing time, labor and costs over traditional methods.

STORAGE CONDITIONS UPON RECEIPT

- Recommended storage for kit components is 2 - 8°C.
- For long-term storage, freeze reconstituted standards and controls at $\leq -20^{\circ}\text{C}$. Avoid multiple (>2) freeze/thaw cycles.
- **DO NOT FREEZE Antibody-Immobilized Beads, Detection Antibody, and Streptavidin-Phycoerythrin.**

REAGENTS SUPPLIED

Note: Store all reagents at 2 – 8 °C

Reagents Supplied	Catalog Number	Volume	Quantity
Pituitary Standard	PIT-8046	lyophilized	1 vial
Pituitary Quality Controls 1 and 2	PIT-6046	lyophilized	2 vials
Set of one 96-Well black Plate with 2 sealers	-----	-----	1 plate 2 sealers
Assay Buffer	L-AB	30 mL	1 bottle
Serum Matrix Note: Contains 0.08% Sodium Azide	PIT-SM	1.5 mL	1 bottle
Bead Diluent	LBD	3.5 mL	1 bottle
10X Wash Buffer Note: Contains 0.05% Proclin	L-WB	30 mL	2 bottles
Pituitary Detection Antibodies	PIT-1046	3.2 mL	1 bottle
Streptavidin-Phycoerythrin	L-SAPE7	5.5 mL	1 bottle
Mixing Bottle	-----	-----	1 bottle

Non-Human Primate Pituitary Hormone Antibody Immobilized Magnetic Beads:

Bead/Analyte Name	Luminex Magnetic Bead Region	Customizable 7 Analytes (20X concentration, 200 µL)	
		Available	Cat. #
ACTH	12	✓	HACTH-MAG
AGRP	15	✓	HAGRP-MAG
CNTF	20	✓	HCNTF-MAG
FSH	26	✓	RFSH-MAG
GH	29	✓	HGH-MAG
LH	33	✓	HLH-MAG
TSH	61	✓	HTSH-MAG

MATERIALS REQUIRED BUT NOT PROVIDED

Reagents

1. Luminex[®] Sheath Fluid (Luminex[®] Catalog #40-50000) or Luminex[®] Drive Fluid (Luminex[®] Catalog #MPXDF-4PK)

Instrumentation / Materials

1. Adjustable Pipettes with Tips capable of delivering 25 μ L to 1000 μ L
2. Multichannel Pipettes capable of delivering 5 μ L to 50 μ L or 25 μ L to 200 μ L
3. Reagent Reservoirs
4. Polypropylene Microfuge Tubes
5. Aluminum Foil
6. Absorbent Pads
7. Rubber Bands
8. Laboratory Vortex Mixer
9. Sonicator (Branson Ultrasonic Cleaner Model # B200 or equivalent)
10. Titer Plate Shaker (Lab-Line Instruments, Model #4625, or equivalent)
11. Luminex[®] 200[™], HTS, FLEXMAP 3D[®], or MAGPIX[®] with xPONENT[®] software by Luminex[®] Corporation
12. Automatic Plate washer for magnetic beads (BioTek[®] ELx405, EMD Millipore Catalog #40-015 or equivalent) or Handheld Magnetic Separation Block (EMD Millipore Catalog #40-285 or equivalent)

Note: If a plate washer or hand held magnetic separation block for magnetic beads is not available, one can use Microtiter filter plate (EMD Millipore) to run the assay with the use of Vacuum Filtration Unit (EMD Millipore Vacuum Manifold Catalog #MSVMHTS00, or equivalent. EMD Millipore Vacuum Pump Catalog #WP6111560 or equivalent.)

SAFETY PRECAUTIONS

- All blood components and biological materials should be handled as potentially hazardous. Follow universal precautions as established by the Centers for Disease Control and Prevention and by the Occupational Safety and Health Administration when handling and disposing of infectious agents.
- Sodium Azide or Proclin has been added to some reagents as a preservative. Although the concentrations are low, Sodium Azide and Proclin may react with lead and copper plumbing to form highly explosive metal azides. Dispose of unused contents and waste in accordance with international, federal, state, and local regulations.

TECHNICAL GUIDELINES

To obtain reliable and reproducible results, the operator should carefully read this entire manual and fully understand all aspects of each assay step before running the assay. The following notes should be reviewed and understood before the assay is set-up.

- FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.
- Do not use beyond the expiration date on the label.
- Do not mix or substitute reagents with those from other lots or sources.

TECHNICAL GUIDELINES (continued)

- The Antibody-Immobilized Beads are light sensitive and must be protected from light at all times. Cover the assay plate containing beads with opaque plate lid or aluminum foil during all incubation steps.
- It is important to allow all reagents to warm to room temperature (20-25°C) before use in the assay.
- Incomplete washing can adversely affect the assay outcome. All washing must be performed with the Wash Buffer provided.
- The standards prepared by serial dilution must be used within 1 hour of preparation. Discard any unused standards except the standard stock which may be stored at $\leq -20^{\circ}\text{C}$ for 1 month and at $\leq -80^{\circ}\text{C}$ for greater than one month.
- If samples fall outside the dynamic range of the assay, further dilute the samples with Assay Buffer or the appropriate diluent and repeat the assay.
- Any unused mixed Antibody-Immobilized Beads may be stored in the Bead Mix bottle at 2-8°C for up to one month.
- During the preparation of the standard curve, make certain to mix the higher concentration well before making the next dilution. Use a new tip with each dilution.
- The plate should be read immediately after the assay is finished. If, however, the plate cannot be read immediately, seal the plate, cover with aluminum foil or an opaque lid, and store the plate at 2-8°C for up to 24 hours. Prior to reading, agitate the plate on the plate shaker at room temperature for 10 minutes. Delay in reading a plate may result in decreased sensitivity for some analytes.
- The titer plate shaker should be set at a speed to provide maximum orbital mixing without splashing of liquid outside the wells. For the recommended plate shaker, this would be a setting of 5-7, which is approximately 500-800 rpm.
- Ensure that the needle probe is clean. This may be achieved by sonication and/or alcohol flushes.
- When reading the assay on Luminex[®] 200™, adjust probe height according to the protocols recommended by Luminex[®] to the kit solid plate using 4 alignment discs. When reading the assay on FLEXMAP 3D[®], adjust probe height according to the protocols recommended by Luminex[®] to the kit solid plate using 1 alignment disc. When reading the assay on MAGPIX[®], adjust probe height according to the protocols recommended by Luminex[®] to the kit solid plate using 2 alignment discs.
- For cell culture supernatants or tissue extraction, use the culture or extraction medium as the matrix solution in blank, standard curve and controls
- For serum / plasma samples, use the matrix provided in the kit in the kit as the matrix solution in blank, standard curve and controls.
- For cell/tissue homogenate, the final cell or tissue homogenate should be prepared in a buffer that has a neutral pH, contains minimal detergents or strong denaturing detergents, and has an ionic strength close to physiological concentration. Avoid debris, lipids, and cell/tissue chunks. Centrifuge samples before use.
- Vortex all reagents well before adding to plate.

SAMPLE COLLECTION AND STORAGE

A. Preparation of Serum Samples:

- Allow the blood to clot for at least 30 minutes before centrifugation for 10 minutes at 1000xg. Remove serum and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$.
- Avoid multiple (>2) freeze/thaw cycles.
- When using frozen samples, it is recommended to thaw the samples completely, mix well by vortexing and centrifuge prior to use in the assay to remove particulates.

B. Preparation of Plasma Samples:

- Plasma collection using EDTA as an anticoagulant is recommended. Centrifuge for 10 minutes at 1000xg within 30 minutes of blood collection. Remove plasma and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$. If measuring ACTH, store plasma samples at -70°C .
- Avoid multiple (>2) freeze/thaw cycles.
- When using frozen samples, it is recommended to thaw the samples completely, mix well by vortexing and centrifuge prior to use in the assay to remove particulates.

C. Preparation of Tissue Culture Supernatant:

- Centrifuge the sample to remove debris and assay immediately or aliquot and store samples at $\leq -20^{\circ}\text{C}$.
- Avoid multiple (>2) freeze/thaw cycles.
- Tissue Culture Supernatant may require a dilution with an appropriate control medium prior to assay.

Note:

- All samples must be stored in polypropylene tubes. **DO NOT STORE SAMPLES IN GLASS.**
- Avoid debris, lipids and cells when using samples with gross hemolysis or lipemia.
- Care must be taken when using heparin as an anticoagulant, since an excess of heparin will provide falsely high values. Use no more than 10 IU heparin per mL of blood collected.

PREPARATION OF REAGENTS FOR IMMUNOASSAY

A. Preparation of Antibody-Immobilized Beads

Sonicate each antibody-bead vial for 30 seconds; vortex for 1 minute. Add 150 μ L from each antibody bead vial to the Mixing Bottle and bring final volume to 3.0 mL with Bead Diluent. Vortex the mixed beads well. Unused portions may be stored at 2-8°C for up to one month. (Note: Due to the composition of magnetic beads, you may notice a slight color in the bead solution. This does not affect the performance of the beads or the kit.)

Example 1: When using 7 antibody-immobilized beads, add 150 μ L from each of the 7 bead sets to the Mixing Bottle. Then add 1.95 mL Bead Diluent.

Example 2: When using 3 antibody-immobilized beads, add 150 μ L from each of the 3 bead sets to the Mixing Bottle. Then add 2.55 mL Bead Diluent.

B. Preparation of Quality Controls

Before use, reconstitute Quality Control 1 and Quality Control 2 with 250 μ L Deionized Water. Invert the vial several times to mix and vortex. Allow the vial to sit for 5-10 minutes. Unused portions may be stored at $\leq -20^{\circ}\text{C}$ for up to one month.

C. Preparation of Wash Buffer

Bring the 10X Wash Buffer to room temperature and mix to bring all salts into solution. Dilute 30 mL of 10X Wash Buffer with 270 mL deionized water. Store unused portions at 2-8°C for up to one month.

D. Preparation of Serum Matrix

This step is required for serum or plasma samples only.

Thaw Serum Matrix completely. Mix well. **Dilute Serum Matrix 1:2 by adding 1.5 mL of Assay Buffer** and mix well. Leftover Serum Matrix should be stored at $\leq -20^{\circ}\text{C}$ for up to one month.

E. Preparation of Pituitary Standard

1.) Prior to use, reconstitute the Pituitary Standard with 250 μ L Deionized Water. Invert the vial several times to mix. Vortex the vial for 10 seconds. Allow the vial to sit for 5-10 minutes. This will be used as Standard 7.


2.) Preparation of Working Standards

Label 6 polypropylene microfuge tubes "Standard 6," "Standard 5," "Standard 4," "Standard 3," "Standard 2," and "Standard 1." Add 150 μ L Assay Buffer to each of the six tubes. Perform 4 times serial dilutions by adding 50 μ L of the "Standard 7" to the "Standard 6" tube, mix well and transfer 50 μ L of the "Standard 6" to the "Standard 5" tube, mix well and transfer 50 μ L of the "Standard 5" to "Standard 4" tube, mix well and transfer 50 μ L of the "Standard 4" to the "Standard 3", mix well and transfer 50 μ L of the "Standard 3" to the "Standard 2" tube, mix well and transfer 50 μ L of the "Standard 2" to the "Standard 1". The 0 Standard (background) will be Assay Buffer.

Preparation of Working Standards:

Reconstituted Standard	Volume of Deionized Water to Add	Volume of Standard to Add
Standard 7	250 μ L	0

Standard	Volume of Assay Buffer to Add	Volume of Standard to Add
Standard 6	150 μ L	50 μ L of Standard 7
Standard 5	150 μ L	50 μ L of Standard 6
Standard 4	150 μ L	50 μ L of Standard 5
Standard 3	150 μ L	50 μ L of Standard 4
Standard 2	150 μ L	50 μ L of Standard 3
Standard 1	150 μ L	50 μ L of Standard 2



After serial dilution, the tubes should have the following concentrations for constructing standard curves.

Standard Tube #	AGRP (pg/mL)	FSH (mIU/mL)	GH (pg/mL)	LH (mIU/mL)	TSH (μ IU/mL)	ACTH (pg/mL)	CNTF (pg/mL)
1	2.4	0.024	2.4	0.049	0.039	3	122
2	10	0.098	10	0.195	0.156	12	488
3	39	0.39	39	0.781	0.625	49	1,953
4	156	1.56	156	3.125	2.5	195	7,813
5	625	6.25	625	12.5	10	781	31,250
6	2,500	25	2,500	50	40	3,125	125,000
7	10,000	100	10,000	200	160	12,500	500,000

IMMUNOASSAY PROCEDURE

- Prior to beginning this assay, it is imperative to read this protocol completely and to thoroughly understand the Technical Guidelines.
- Allow all reagents to warm to room temperature (20-25°C) before use in the assay.
- Diagram the placement of Standards, 0 (Background), Std 1, Std 2, Std 3, Std 4, Std 5, Std 6, and Std 7, Controls 1 and 2, and samples on Well Map Worksheet in a vertical configuration. (Note: Most instruments will only read the 96-well plate vertically by default.) It is recommended to run the assay in duplicate.
- If using a filter plate, set on a plate holder at all times during reagent dispensing and incubation steps so that the bottom of the plate does not touch any surface.

1. Add 200 µL of Assay Buffer into each well of the plate. Seal and mix on a plate shaker for 10 minutes at room temperature (20-25°C).
2. Decant Assay Buffer and remove the residual amount from all wells by inverting the plate and tapping it smartly onto absorbent towels several times.
3. Add 25 µL of 1:2 diluted Matrix Solution (when measuring serum or plasma samples) or appropriate culture media (when measuring culture samples) in Background, Standards, and Quality Control wells.
4. Add 25 µL Assay Buffer in Sample wells
5. Add 25 µL of Assay Buffer to the Background wells. Add 25 µL of each Standard or Control into the appropriate wells.
6. Add 25 µL of samples to the Sample wells.
7. Vortex Bead Bottle and add 25 µL of the prepared Beads to each well. (Note: during addition of the Beads, shake beads intermittently to avoid settling)
8. Seal the plate with a plate sealer (wrap the plate with foil if not using foil pate sealer) and incubate with agitation on a plate shaker for overnight incubation at 4°C (16-18 hr).

Add 200 µL Assay Buffer
per well



Shake 10 min, RT

Decant

- Add 25 µL Matrix Solution (or appropriate media) to Background, Standard and Quality Controls.
- Add 25 µL Assay Buffer to Sample wells.
- Add 25µL Assay Buffer to background wells and 25µL Standard, and 25µL Controls to Standard and Control wells, respectively.
- Add 25µL samples to Sample Wells.
- Add 25 µL Beads to each well



Incubate overnight
at 4°C with agitation
on a plate shaker

9. Gently remove well contents and wash plate 2 times following instructions listed in the **PLATE WASHING** section.
10. Add 25 μL of Detection Antibodies into each well. (Note: allow the Detection Antibodies to warm to room temperature prior to addition.)
11. Seal, cover with foil, and incubate with agitation on a plate shaker for 60 minutes at room temperature (20-25°C) **DO NOT ASPIRATE AFTER INCUBATION.**
12. Add 25 μL Streptavidin-Phycoerythrin to each well containing the 25 μL of Detection Antibodies.
13. Seal, cover with foil and incubate with agitation on a plate shaker for 30 minutes at room temperature (20-25°C).
14. Gently remove well contents and wash plate 2 times following instructions listed in the **PLATE WASHING** section.
15. Add 100 μL of Sheath Fluid (or Drive Fluid if using MAGPIX[®]) to all wells. Resuspend the beads on a plate shaker for 5 minutes.
16. Run plate on Luminex[®] 200[™], HTS, FLEXMAP 3D[®] or MAGPIX[®] with xPONENT[®] software.
17. Save and analyze the median Fluorescent Intensity (MFI) data using a weighted 5-parameter logistic or spline curve-fitting method for calculating analyte concentrations in samples.



Remove well contents. Wash 2X with 200 μL Wash Buffer

Add 25 μL Detection Antibody per well



Incubate 60 min at RT

Do Not Aspirate

Add 25 μL Streptavidin-Phycoerythrin per well

Incubate for 30 minute at RT



Remove well contents and Wash 2X with 200 μL Wash Buffer

Add 100 μL Sheath Fluid or Drive Fluid per well

Read on Luminex[®] (100 μL , 50 Beads per Bead set)

PLATE WASHING

1.) Solid Plate

If using a solid plate, use either a hand-held magnet or magnetic plate washer.

A.) For hand-held magnet, rest plate on magnet for 60 seconds to allow complete settling of magnetic beads. Remove well contents by gently decanting the plate in an appropriate waste receptacle and gently tapping on absorbent pads to remove residual liquid. Wash plate with 200 μ L of Wash Buffer by removing plate from magnet, adding Wash Buffer, shaking for 30 seconds, reattaching to magnet, letting beads settle for 60 seconds and removing well contents as previously described after each wash. Repeat wash steps as recommended in Assay Procedure.

B.) For magnetic plate washer, let plate "soak" on magnet for 60 seconds to allow complete settling of the magnetic beads. Remove well contents by aspiration. Wash plate with 200 μ L/well of Wash Buffer, letting beads "soak" for 60 seconds and removing Wash Buffer by aspiration after each wash. Repeat wash steps as recommended in Assay Procedure. **Note:** If using the recommended plate washer for magnetic beads (Bio-Tek ELx405) follow the appropriate equipment settings outlined in **EQUIPMENT SETTINGS**.

2.) Filter Plate (EMD Millipore Cat# MX-PLATE)

If using a filter plate, use a vacuum filtration manifold to remove well contents. Wash plate 200 μ L/well of Wash Buffer, removing Wash Buffer by vacuum filtration after each wash. Repeat wash steps as recommended in the Assay Procedure.

EQUIPMENT SETTINGS

Bio-Tek® ELx405:

The general recommended wash protocol (Link Protocol) is as follows:

Soak Program: Wash Program:

Soak → Aspirate→Dispense→Soak→Aspirate→Dispense→Soak→Aspirate

1.) Soak program:

1. Soak duration: 60 sec
2. Shake before soak?: NO

2.) Wash program:

Method:

1. Number of cycles: 2
2. soak/shake: YES
3. Soak duration: 60 sec
4. Shake before soak: NO
5. Prime after soak: NO

Dispense:

1. Dispense volume: 200 µL/well
2. Dispense flow rate: 5
3. Dispense height: 130 (16.51 mm)
4. Horizontal disp pos: 00 (0 mm)
5. Bottom Wash first?: NO
6. Prime before start?: NO

Aspiration:

1. Aspirate height: 35 (4.445 mm)
2. Horizontal Asp Pos: 30 (1.372 mm)
3. Aspiration rate: 06 (15.0 mm/sec)
4. Aspiration delay: 0
5. Crosswise Aspir: NO
6. Final Aspir: YES
7. Final Aspir delay: 0 (0 msec)

3.) Link program: (**Note:** this is the program to use during actual plate washing).
Link together the Soak and Wash programs outlined above.

Note: After the final aspiration, there will be approximately 25 µL of residual Wash Buffer in each well. This is expected when using the BioTek Plate washer and this volume does not need to be aspirated from the plate.

If using an automatic plate washer other than BioTek ELx405, please refer to the manufacturer's recommendations for programming instructions.

EQUIPMENT SETTINGS (continued)

Luminex[®] 200[™], HTS, FLEXMAP 3D[®] and MAGPIX[®] with xPONENT[®] software:

These specifications are for the Luminex[®] 200[™], Luminex[®] HTS, Luminex[®] FLEXMAP 3D[®] and Luminex[®] MAGPIX[®] with xPonent[®] software. Luminex[®] instruments with other software (e.g. MasterPlex[®], STarStation, LiquiChip, Bio-Plex Manager[™], LABScan[™] 100) would need to follow instrument instructions for gate settings and additional specifications from the vendors for reading Luminex[®] Magnetic Beads.

For magnetic bead assays, the Luminex[®] 200[™] and HTS instruments must be calibrated with the xPONENT[®] 3.1 compatible Calibration Kit (EMD Millipore Catalog #40-275) and performance verified with the Performance Verification Kit (EMD Millipore Catalog #40-276). The Luminex[®] FLEXMAP 3D[®] instrument must be calibrated with the FLEXMAP 3D[®] Calibrator Kit (EMD Millipore Catalog #40-028) and performance verified with the FLEXMAP 3D[®] Performance Verification Kit (EMD Millipore Catalog #40-029). The Luminex[®] MAGPIX[®] instrument must be calibrated with the MAGPIX[®] Calibration Kit (EMD Millipore Catalog #40-049) and performance verified with the MAGPIX[®] Performance Verification Kit (EMD Millipore Catalog #40-050).

NOTE: These assays cannot be run on any instruments using Luminex[®] IS 2.3 or Luminex[®] 1.7 software.

The Luminex[®] probe height must be adjusted to the plate provided in the kit. Please use Cat# MAG-PLATE, if additional plates are required for this purpose.

Events:	50, per bead	
Sample Size:	100 μ L	
Gate Settings:	8,000 to 15,000	
Reporter Gain:	Default (low PMT)	
Time Out:	60 seconds	
Bead Set:	7-Plex Beads	
	ACTH Beads	12
	AGRP Beads	15
	CNTF Bead	20
	FSH Bead	26
	GH Bead	29
	LH Bead	33
	TSH Bead	61

QUALITY CONTROLS

The ranges for each analyte in Quality Control 1 and 2 are provided on the card insert or can be located at the EMD MILLIPORE website www.millipore.com/techlibrary/index.do using the Catalog number as the keyword.

ASSAY CHARACTERISTICS

Assay Sensitivities (minimum detectable concentrations)

Minimum Detectable Concentration (MinDC) defines the reportable range of the assay. It is a measure of the true limits of detection for an assay and is mathematically determined.

Analyte	Overnight Protocol (N = 6 assays)	
	Mean MinDC	Mean MinDc + 2SD
ACTH (pg/mL)	0.79	1.42
AGRP (pg/mL)	0.80	1.81
CNTF (pg/mL)	37.42	100.94
FSH (mIU/mL)	0.012	0.02
GH (pg/mL)	2.03	4.03
LH (mIU/mL)	0.013	0.02
TSH (μ IU/mL)	0.005	0.01

Precision

Intra-assay precision is generated from the mean of the %CV's from 16 reportable results across two different concentration of analytes in one experiment. Inter-assay precision is generated from the mean of the %CV's from two reportable results each for two different concentrations of cytokine different experiments.

Analyte	Intra-Assay CV %	Inter-Assay CV %
ACTH	<10	<15
AGRP	<10	<15
CNTF	<10	<15
FSH	<10	<15
GH	<10	<15
LH	<10	<15
TSH	<10	<15

Accuracy

Spike Recovery: The data represents mean recovery of 3 levels of spiked standards using 5 independent serum samples.

Analyte	Spike and Recovery %
ACTH	112
AGRP	93
CNTF	71
FSH	98
GH	106
LH	70
TSH	113

Cross-Reactivity

The antibody pairs in the panel are specific only to the desired analyte and exhibit no or negligible (<2%) cross-reactivity with other analytes in the panel.

TROUBLESHOOTING GUIDE

Problem	Probable Cause	Solution
Insufficient Bead Count	<p>Plate Washer aspirate height set too low</p> <p>Bead mix prepared inappropriately</p> <p>Samples cause interference due to particulate matter or viscosity</p> <p>Probe height not adjusted correctly</p>	<p>Adjust aspiration height according to manufacturers' instructions.</p> <p>Sonicate bead vials and vortex just prior to adding to bead mix bottle according to protocol. Agitate bead mix intermittently in reservoir while pipetting this into the plate.</p> <p>See above. Also sample probe may need to be cleaned with Alcohol flush, Back flush and washes; or if needed probe should be removed and sonicated.</p> <p>When reading the assay on Luminex[®] 200™, adjust probe height according to the protocols recommended by Luminex[®] to the kit solid plate using 4 alignment discs. When reading the assay on FLEXMAP 3D[®], adjust probe height according to the protocols recommended by Luminex[®] to the kit solid plate using 1 alignment disc. When reading the assay on MAGPIX[®], adjust probe height according to the protocols recommended by Luminex[®] to the kit solid plate using 2 alignment discs.</p>
Background is too high	<p>Background wells were contaminated</p> <p>Matrix used has endogenous analyte or interference</p> <p>Insufficient washes</p>	<p>Avoid cross-well contamination by using sealer appropriately, and pipeting with Multichannel pipets without touching reagent in plate.</p> <p>Check matrix ingredients for cross reacting components (e.g. interleukin modified tissue culture medium).</p> <p>Increase number of washes.</p>
Beads not in region or gate	<p>Luminex[®] not calibrated correctly or recently</p> <p>Gate Settings not adjusted correctly</p> <p>Wrong bead regions in protocol template</p> <p>Incorrect sample type used</p> <p>Instrument not washed or primed</p>	<p>Calibrate Luminex[®] based on Instrument Manufacturer's instructions, at least once a week or if temperature has changed by >3°C.</p> <p>Some Luminex[®] instruments (e.g. Bioplex) require different gate settings than those described in the Kit protocol. Use Instrument default settings.</p> <p>Check kit protocol for correct bead regions or analyte selection.</p> <p>Samples containing organic solvents or if highly viscous should be diluted or dialyzed as required.</p> <p>Prime the Luminex[®] 4 times to rid of air bubbles, wash 4 times with sheath fluid or water if there is any remnant alcohol or sanitizing liquid.</p>

Problem	Probable Cause	Solution
Beads not in region or gate (continued)	Beads were exposed to light	Keep plate and bead mix covered with dark lid or aluminum foil during all incubation steps.
Signal for whole plate is same as background	Incorrect or no Detection Antibody was added Streptavidin-Phycoerythrin was not added	Add appropriate Detection Antibody and continue. Add Streptavidin-Phycoerythrin according to protocol. If Detection Antibody has already been removed, sensitivity may be low.
Low signal for standard curve	Detection Antibody may have been removed prior to adding Streptavidin Phycoerythrin Incubations done at inappropriate temperatures, timings or agitation	May need to repeat assay if desired sensitivity not achieved. Assay conditions need to be checked.
Signals too high, standard curves are saturated	Calibration target value set too high Plate incubation was too long with standard curve and samples	With some Luminex [®] Instrument (e.g. Bio-plex) Default target setting for RP1 calibrator is set at High PMT. Use low target value for calibration and reanalyze plate. Use shorter incubation time.
Sample readings are out of range	Samples contain no or below detectable levels of analyte Samples contain analyte concentrations higher than highest standard point. Standard curve was saturated at higher end of curve.	If below detectable levels, it may be possible to use higher sample volume. Check with tech support for appropriate protocol modifications. Samples may require dilution and reanalysis for just that particular analyte. See above.
High Variation in samples and/or standards	Multichannel pipet may not be calibrated Plate washing was not uniform Samples may have high particulate matter or other interfering substances Plate agitation was insufficient Cross well contamination	Calibrate pipets. Confirm all reagents are removed completely in all wash steps. See above. Plate should be agitated during all incubation steps using a vertical plate shaker at a speed where beads are in constant motion without causing splashing. Check when reusing plate sealer that no reagent has touched sealer.

Problem	Probable Cause	Solution
High Variation in samples and/or standards	Cross well contamination	Care should be taken when using same pipet tips that are used for reagent additions and that pipet tip does not touch reagent in plate.
FOR FILTER PLATES ONLY		
Filter plate will not vacuum	<p>Vacuum pressure is insufficient</p> <p>Samples have insoluble particles</p> <p>High lipid concentration</p>	<p>Increase vacuum pressure such that 0.2mL buffer can be suctioned in 3-5 seconds.</p> <p>Centrifuge samples just prior to assay setup and use supernatant.</p> <p>After centrifugation, remove lipid layer and use supernatant.</p>
Plate leaked	<p>Vacuum Pressure too high</p> <p>Plate set directly on table or absorbent towels during incubations or reagent additions</p> <p>Insufficient blotting of filter plate bottom causing wicking</p> <p>Pipette touching plate filter during additions</p> <p>Probe height not adjusted correctly</p> <p>Sample too viscous</p>	<p>Adjust vacuum pressure such that 0.2mL buffer can be suctioned in 3-5 seconds. May need to transfer contents to a new (blocked) plate and continue.</p> <p>Set plate on plate holder or raised edge so bottom of filter is not touching any surface.</p> <p>Blot the bottom of the filter plate well with absorbent towels after each wash step.</p> <p>Pipette to the side of plate.</p> <p>Adjust probe to 3 alignment discs in well H6.</p> <p>May need to dilute sample.</p>

REPLACEMENT REAGENTS

	Cat #
Pituitary Standard	PIT-8046
Pituitary Quality Controls	PIT-6046
Pituitary Detection Antibodies	PIT-1046
Serum Matrix	PIT-SM
Bead Diluent	LBD
Assay Buffer	L-AB
Streptavidin-Phycoerythrin	L-SAPE7
Set of two 96-Well Black plates with 4 sealers	MAG-PLATE
10X Wash Buffer	L-WB

Antibody-Immobilized Beads

Analyte	Bead #	Cat. #
ACTH	12	HACTH-MAG
AGRP	15	HAGRP-MAG
CNTF	20	HCNTF-MAG
FSH	26	RFSH-MAG
GH	29	HGH-MAG
LH	33	HLH-MAG
TSH	61	HTSH-MAG

ORDERING INFORMATION

To place an order:

To assure the clarity of your custom kit order, please FAX the following information to our customer service department:

Include:

- Your name, telephone and/or fax number
- Customer account number
- Shipping and billing address
- Purchase order number
- Catalog number and description of product
- Quantity of kits
- Selection of MILLIPLEX[®] Analytes

FAX: (636) 441-8050

Toll-Free US: (866) 441-8400
(636) 441-8400

Mail Orders: EMD Millipore Corporation
6 Research Park Drive
St. Charles, Missouri 63304 U.S.A.

For International Customers:

To best serve our international customers in placing an order or obtaining additional information about MILLIPLEX[®] MAP products, please contact your multiplex specialist or sales representative or email our European Customer Service at customerserviceEU@merckgroup.com.

Conditions of Sale

For Research Use Only. Not for Use in Diagnostic Procedures.

Material Safety Data Sheets (MSDS)

Material Safety Data Sheets for EMD Millipore products may be ordered by fax or phone or through our website at www.emdmillipore.com/techlibrary/index.do.

Technical Services

<http://www.emdmillipore.com/techservices>

To contact by phone

For North America: Toll-Free US: 1-(800) 221-1975 or 1-(781) 533-8045

Outside North America, contact your local office

<http://www.emdmillipore.com/offices>

WELL MAP

	1	2	3	4	5	6	7	8	9	10	11	12
A	0 Standard (Background)	Standard 4	QC-1 Control	Etc.								
B	0 Standard (Background)	Standard 4	QC-1 Control									
C	Standard 1	Standard 5	QC-2 Control									
D	Standard 1	Standard 5	QC-2 Control									
E	Standard 2	Standard 6	Sample 1									
F	Standard 2	Standard 6	Sample 1									
G	Standard 3	Standard 7	Sample 2									
H	Standard 3	Standard 7	Sample 2									