

DNA Isolation from **1 ml** or **2 ml** plasma/serum. (50 isolations)

- 1) Add **30 µl Plasma Condition Solution** to **1 ml** plasma or **60 µl** to **2 ml** of plasma and vortex for 5 seconds.
- 2) Add **140 µl** of **SubX[®]** solution to **1 ml** of human plasma or **280 µl** of **SubX[™]** solution to **2 ml** of human plasma, close the tube and immediately vortex for 10-15 seconds. Turbidity occurs. Incubate 3 min at room temperature; vortex 3-4 times during incubation.
- 3) Vortex **Binding Matrix** for 15 sec; open the tube and add **20 µl** of **Binding Matrix** slurry to **SubX[™]**-plasma mixture, close the tube and incubate 3 min at room temperature; briefly (~10 sec) vortex 3-4 times during incubation.
- 4) Centrifuge the tube for 1 min at 14,000 rpm, aspirate the supernatant carefully by pipetting taking care not to disturb the **Binding Matrix**, discard supernatant.
- 5) Add **0.5 ml Wash-1** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 6) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 7) Add **0.5 ml Wash-G** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 8) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 9) Add **0.6 ml Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 10) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 11) Add **0.6 ml Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 12) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 13) **ESSENTIAL step!** Briefly centrifuge the tube for 15 sec at 14,000 rpm to remove residual traces of washing buffer from the tube wells, then aspirate traces of Wash-2 solution by using a 10-20 µl micropipette tip and discard.
- 14) Add 30-50 µl of **Elution** solution to the pellet, tightly close the tube and incubate at +60°C for 5 min. Briefly (~10 sec) vortex the tube 5-6 times during incubation.
- 15) Centrifuge the tube for 1 min at 14,000 rpm, carefully aspirate the cfDNA containing supernatant, transfer to a new tube and **save**.
- 16) **(Optional)** Repeated elution can yield additional 10-20% of DNA.
- 17) Keep the cfDNA containing supernatant at +4°C for next-day use or at -20°C for long-term storage.

DNA recovery: Our kit allows for isolation of DNA fragments under 50 bp and even less. Typical DNA recovery ranges from 1 to 100 ng/ml of plasma. The yield varies depending on the sample source and donor health conditions.

Contact Us

Capital Biosciences, Inc.
900 Clopper Rd, Suite 120
Gaithersburg, MD 20878
Phone: **1-800-475-2812**

Email: info@capitalbiosciences.com

Web: www.capitalbiosciences.com

SubX[™] Plasma cfDNA isolation kit

INSTRUCTION MANUAL (v1.5)

Catalog No. CFDNA-0050
For Research Use Only

Product Description

SubX[™] kit is designed for isolation of circulating cell-free DNA (cfDNA) directly from liquid biopsies (plasma or serum) without Proteinase K. Our technology is based on the use of proprietary bi-functional substance (SubX[™]) that binds DNA under physiological conditions (e.g. directly in biological liquids), followed by adsorption of [DNA-SubX[™]] complex on a solid phase matrix. Since SubX[™] captures DNA via phosphate residues (groups) it allows for elimination of bias related to both AT/CG content and DNA fragments length, thus improving extraction efficacy and accuracy of downstream applications. All currently available commercial kits are [silica-chaotropic salt]-based and exhibit bias for either short or long DNA fragments, as well as for the GC content. In addition, such systems suffer from lot-to-lot minor surface differences of silica that results in variations of extraction efficiency of small DNA fragments. Our approach eliminates these problems. Specially designed solid-phase matrix does not require addition of (unwanted) exogenous RNA to prevent non-specific DNA adsorption. Universal conditioning buffer makes our kit compatible with all tested commercial serum/plasma storage systems/solutions

Kit components, shipping and storage conditions:

- **SubX[™] Solution** (2x7 ml). Ready-to-use. **Store refrigerated +4°C.**
- **Plasma Condition Solution** (1.5 ml) Ready-to-use. Store at room temperature.
- **Binding Matrix** (1 ml). Ready-to-use. Store at room temperature in a dark place.
- **Wash-1 Solution** (30 ml). Ready-to-use. Store at room temperature in a dark place.
- **Wash-G Solution** (12 ml). **Add 18 ml of 96% molecular grade ethanol before use.**
- **Wash-2 Solution** (2x11 ml). **Add 55 ml of 96% molecular grade ethanol before use.** Store at room temperature.
- **Elution Buffer:** (5 ml); Ready-to-use. **Store refrigerated +4°C.**

Kits are shipped at ambient temperatures.

Important Notes (Read this first!!!)

- To prevent contamination of circulating cell free DNA with cellular DNA centrifuge plasma for at least 1-2 min at 14000 rpm before use, carefully aspirate supernatant not to disturb possible cell debris pellet and transfer into new tube.
- Optimal results are achieved when Wash-2 solution is diluted with molecular grade 96% ethanol.
- We recommend using Posi-Click 1.5 ml Eppendorf tubes (Denville Scientific, Inc., Cat. C2170) to prevent material loss during vortexing due to possible cap opening.
- Binding Matrix can be separated from the liquid phase with magnetic stand such as DynaMag[™]-2 Magnet. **Very important: always pellet Binding Matrix by centrifugation after step 3 of the protocol.** If using magnetic separation protocol, centrifugation for a short time is necessary to collect all materials from the tube cap.

DNA Isolation from **0.2 ml** plasma/serum (100 isolations)

- 1) Add **6 µl Plasma Condition Solution** and vortex for 5 seconds.
- 2) Add **30 µl** of **SubX™** solution to **0.2 ml** of human plasma, close the tube and immediately vortex for 10-15 seconds. Turbidity occurs. Incubate 3 min at room temperature; vortex 3-4 times during incubation.
- 3) Vortex **Binding Matrix** for 15 sec; open the tube and add **10 µl** of **Binding Matrix** slurry to **SubX™**-plasma mixture, close the tube and incubate 3 min at room temperature; briefly (~10 sec) vortex 3-4 times during incubation.
- 4) Centrifuge the tube for 1 min at 14,000 rpm, aspirate the supernatant carefully by pipetting taking care not to disturb the **Binding Matrix**, discard supernatant.
- 5) Add **0.3 ml Wash-1** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 6) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 7) Add **0.3 ml Wash-G** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 8) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 9) Add **0.4 ml Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 10) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 11) Add **0.4 ml Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 12) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 13) **ESSENTIAL step!** Briefly centrifuge the tube for 15 sec at 14,000 rpm to remove residual traces of washing buffer from the tube wells, then aspirate traces of Wash-2 solution by using a 10-20 µl micropipette tip and discard.
- 14) Add 20-50 µl of **Elution** solution to the pellet, tightly close the tube and incubate at +60°C for 5 min. Briefly (~10 sec) vortex the tube 5-6 times during incubation.
- 15) Centrifuge the tube for 1 min at 14,000 rpm, carefully aspirate the cfDNA containing supernatant, transfer to a new tube and **save**.
- 16) **(Optional)** Repeated elution can yield additional 10-20% of DNA.
- 17) Keep the cfDNA containing supernatant at +4°C for next-day use or at -20°C for long-term storage.

DNA Isolation from **0.5 ml** plasma/serum (65 isolations)

- 1) Add **15 µl Plasma Condition Solution** and vortex for 5 seconds.
- 2) Add **70 µl** of **SubX™** solution to **0.5 ml** of human plasma, close the tube and immediately vortex for 10-15 seconds. Turbidity occurs. Incubate 3 min at room temperature; vortex 3-4 times during incubation.
- 3) Vortex **Binding Matrix** for 15 sec; open the tube and add **15 µl** of **Binding Matrix** slurry to **SubX™**-plasma mixture, close the tube and incubate 3 min at room temperature; briefly (~10 sec) vortex 3-4 times during incubation.
- 4) Centrifuge the tube for 1 min at 14,000 rpm, aspirate the supernatant carefully by pipetting taking care not to disturb the **Binding Matrix**, discard supernatant.
- 5) Add **0.4 ml Wash-1** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 6) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 7) Add **0.4 ml Wash-G** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 8) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 9) Add **0.5 ml Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 10) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 11) Add **0.5 ml Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 12) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 13) **ESSENTIAL step!** Briefly centrifuge the tube for 15 sec at 14,000 rpm to remove residual traces of washing buffer from the tube wells, then aspirate traces of Wash-2 solution by using a 10-20 µl micropipette tip and discard.
- 14) Add 20-50 µl of **Elution** solution to the pellet, tightly close the tube and incubate at +60°C for 5 min. Briefly (~10 sec) vortex the tube 5-6 times during incubation.
- 15) Centrifuge the tube for 1 min at 14,000 rpm, carefully aspirate the cfDNA containing supernatant, transfer to a new tube and **save**.
- 16) **(Optional)** Repeated elution can yield additional 10-20% of DNA.
- 17) Keep the cfDNA containing supernatant at +4°C for next-day use or at -20°C for long-term storage.