

Exosome and DNA Isolation from **1 ml** or **2 ml** plasma/serum. (up to 50 isolations)

- 1) Add **30 µl** **Stabilization (Plasma Condition) Solution** to **1 ml** plasma or **60 µl** to **2 ml** of plasma and vortex for 5 seconds.
- 2) Add **200 µl** of **SubX™** solution to **1 ml** of human plasma or **400 µl** of **SubX™** solution to **2 ml** of human plasma, close the tube and immediately vortex for 10-15 seconds. Turbidity may occur. Incubate 3-5 min at room temperature; vortex 3-4 times during incubation.
- 3) Vortex **Binding Matrix (BM)** for 15 sec; open the tube and add **20 µl** of **Binding Matrix** slurry to **SubX™**-plasma mixture, close the tube and incubate 3-5 min at room temperature; briefly (~10 sec) vortex 3-4 times during incubation.
- 4) Centrifuge the tube for 5 min at 14,000 rpm, aspirate the supernatant carefully by pipetting taking care not to disturb the **Binding Matrix**, discard supernatant.
- 5) Discard supernatant, invert tube onto paper towel and let remaining liquid to drain or aspirate carefully by pipette
- 6) Wash [BM-SubX] pellet with 200 µl of **Wash-Ex** buffer
- 7) Centrifuge the tube for 3 min at 14,000 rpm and discard supernatant
- 8) Resuspend pellet in **0.1-0.2 ml** **ERB (Exosome Reconstitution Buffer)**
- 9) Centrifuge at 14000 rpm for 1-5 min to pellet [SubX-DNA-BM]
- 10) Transfer Exosome supernatant into new tube and freeze
- 11) Resuspend [SubX-DNA-BM] pellet in **0.5 ml** **Wash-1** solution, mix and incubate 3-4 min
- 12) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 13) Add **0.5 ml** **Wash-G** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 14) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 15) Add **0.6 ml** **Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 16) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 17) Add **0.6 ml** **Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex 3-4 times for 3 min;
- 18) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 19) **ESSENTIAL step!** Briefly centrifuge the tube for 15 sec at 14,000 rpm to remove residual traces of washing buffer from the tube wells, then aspirate traces of Wash-2 solution by using a 10-20 µl micropipette tip and discard.
- 20) Add 30-50 µl of **DNA Elution Buffer** to the pellet, tightly close the tube and incubate at +60°C for 5 min. Briefly (~10 sec) vortex the tube 5-6 times during incubation.
- 21) Centrifuge the tube for 1 min at 14,000 rpm, carefully aspirate the cfDNA containing supernatant, transfer to a new tube and **save**.
- 22) **(Optional)** Repeated elution can yield additional 10-20% of DNA.
- 23) Keep the cfDNA containing supernatant at +4°C for next-day use or at -20°C for long-term storage.

DNA recovery: Our kit allows for isolation of DNA fragments under 50 bp and even less. Typical DNA recovery ranges from 1 to 100 ng/ml of plasma. The yield varies depending on the sample source and donor health conditions.

SubX™ Exo-cfDNA isolation kit

INSTRUCTION MANUAL (v1.5)

Catalog No. EXOD-0050

For Research Use Only

Product Description

SubX™ Exo-cfDNA kit is designed for isolation of exosomes and circulating cell-free DNA (cfDNA) directly from liquid biopsies (plasma or serum) without Proteinase K. Our technology is based on the use of proprietary bi-functional substance (SubX™) that binds DNA under physiological conditions (e.g. directly in biological liquids), followed by adsorption of [DNA- SubX™] complex on a solid phase matrix. The unique property of SubX™ is that both ends of the molecule bear lipid binding groups. This feature allows each molecule of SubX to anchor two exosomes (i.e. dimerize). Excess of SubX™ in the solution results in aggregation of up to 10-15 exosomes. Since SubX™ captures DNA via phosphate residues (groups) it allows for elimination of bias related to both AT/CG content and DNA fragments length, thus improving extraction efficacy and accuracy of downstream applications. All currently available commercial kits are [silica-chaotropic salt]-based and exhibit bias for either short or long DNA fragments, as well as for the GC content. In addition, such systems suffer from lot-to-lot minor surface differences of silica that results in variations of extraction efficiency of small DNA fragments. Our approach eliminates these problems. Specially designed solid-phase matrix does not require addition of exogenous RNA to prevent non-specific DNA adsorption. Universal conditioning buffer makes our kit compatible with all tested commercial serum/plasma storage systems/solutions

Kit components, shipping and storage conditions:

Plasma Condition Solution (2x 1.5 ml) Ready-to-use. Store at room temperature.

Urine Condition Solution (0.5 ml) Ready-to-use. Store at room temperature.

SubX™ Solution (2x7 ml). Ready-to-use. Store at room temperature

Wash-Ex Solution (10 ml) Ready-to-use. Store at room temperature in a dark place.

Exosome Reconstitution Buffer (ERB) - 5 ml. Store at room temperature.

Binding Matrix (1 ml). Ready-to-use. Store at room temperature in a dark place.

Wash-1 Solution (30 ml). Ready-to-use. Store at room temperature in a dark place.

Wash-G Solution (10 ml). **Add 15 ml of 96% molecular grade ethanol before use.**

Wash-2 Solution (2x11 ml). **Add 55 ml of 96% molecular grade ethanol before use.**

Elution Buffer: (2x1.5 ml); Ready-to-use. Store at room temperature.

Important Notes (Read this first!!!)

- To prevent contamination of circulating cell free DNA with cellular DNA centrifuge plasma for at least 1-2 min at 14000 rpm before use, carefully aspirate supernatant not to disturb possible cell debris pellet and transfer into new tube.
- If working with plasma from **EDTA blood** collection tubes, e.g., Streck Cell-Free DNA BCT™, it is important to add **Plasma Condition Solution** (1/30 volume) before adding SubX™ Solution. One can skip step 1 if working with plasma from **citrated blood** collection.
- Optimal results are achieved when Wash-G and Wash-2 solution is diluted with molecular grade 96% ethanol.
- We recommend using Posi-Click 1.5 ml Eppendorf tubes (Denville Scientific, Inc., Cat. C2170) to prevent material loss during vortexing due to possible cap opening.
- Binding Matrix can be separated from the liquid phase with magnetic stand such as DynaMag™-2 Magnet. **Very important: always pellet Binding Matrix by centrifugation after step 3 of the protocol.** If using magnetic separation protocol, centrifugation for a short time is necessary to collect all materials from the tube cap.

DNA and exosome isolation protocol from **large volume** of urine:

- 1) Collect urine in 50 ml falcon tubes
- 2) Add 10 μ l **Urine Condition Solution per 5 ml of urine**
- 3) Centrifuge 15 min at 3500 rpm (~2400 g) to pellet cell debris
- 4) Transfer supernatant into fresh 40 ml Nalgene centrifugation tubes (e.g., 2x 25 ml)
- 5) Add **1 ml SubX™** per 5 ml supernatant, mix and incubate 30-40 min at room temperature
- 6) Add **BM (30 μ l)**, mix and incubate 3-4 min
- 7) Centrifuge 10 min at 1000 rpm (~14000 g) in high-speed centrifuge
- 8) Discard supernatant, invert tube onto paper towel and let remaining liquid to drain
- 9) Wash **[BM-SubX]** pellet with 300 μ l of **Wash-Ex** buffer. Resuspend the pellet by pipetting and transfer into 1.5 ml Eppendorf tube.
- 10) Centrifuge the tube for 5 min at 14,000 rpm and discard supernatant
- 11) Resuspend pellet in **0.2-0.5 ml ERB (Exosome Reconstitution Buffer)**
- 12) Centrifuge at 14000 rpm for 5 min to pellet **[SubX-DNA-BM]**
- 13) Aspirate Exosome supernatant into new tube and freeze at -80C
- 14) Resuspend **[SubX-DNA-BM]** pellet in **0.5 ml Wash-1** solution, mix and incubate 3-4 min
- 15) Centrifuge for 1 min at 14000 rpm
- 16) Discard supernatant
- 17) Resuspend **[SubX-DNA-BM]** pellet in **0.5 ml Wash-G** solution, mix and incubate 3-4 min
- 18) Centrifuge for 1 min at 14000 rpm
- 19) Discard supernatant
- 20) Resuspend pellet in **0.8 ml Wash-2**, mix and incubate 3-4 min
- 21) Centrifuge for 1 min at 14000 rpm
- 22) Discard supernatant
- 23) Resuspend pellet in **0.8 ml Wash-2**, mix and incubate 3-4 min
- 24) Centrifuge for 1 min at 14000 rpm
- 25) Discard supernatant
- 26) Centrifuge 15 sec at 14000 rpm
- 27) Aspirate leftover using 10-20 μ l tip
- 28) Resuspend pellet in **50 – 200 μ l DNA Elution Buffer**
- 29) Elution: Incubate at 60°C for 5 minutes at 1200 rpm on Thermomixer
- 30) Centrifuge at 14000 rpm for 1 min
- 31) Aspirate DNA supernatant and transfer to new tube and fr

DNA Isolation from exosomes.

1. Dilute **0.1 ml** exosomes with **0.9 ml** PBS.
2. Add **140 μ l** of **SubX™** solution, close the tube and immediately vortex for 10-15 seconds.
3. Incubate **30 min** at room temperature; vortex 3-4 times during incubation.
4. Add **BM (10 μ l)**, mix and incubate 3-4 min
5. Centrifuge at 10,000g for 5 min
6. Discard supernatant
7. Resuspend **[SubX-DNA-BM]** pellet in **0.3 ml Wash-1** solution, mix and incubate 3-4 min
8. Centrifuge for 1 min at 14000 rpm
9. Discard supernatant
10. Resuspend **[SubX-DNA-BM]** pellet in **0.3 ml Wash-G** solution, mix and incubate 3-4 min
11. Centrifuge for 1 min at 14000 rpm
12. Discard supernatant
13. Resuspend pellet in **0.6 ml Wash-2**, mix and incubate 3-4 min
14. Centrifuge for 1 min at 14000 rpm
15. Discard supernatant
16. Resuspend pellet in **0.6 ml Wash-2**, mix and incubate 3-4 min
17. Centrifuge for 1 min at 14000 rpm
18. Discard supernatant
19. Centrifuge 15 sec at 14000 rpm
20. Aspirate leftover using 10-20 μ l tip
21. Resuspend pellet in **50 – 200 μ l DNA Elution Buffer**
22. Elution: Incubate at 60°C for 5 minutes at 1200 rpm on Thermomixer
23. Centrifuge at 14000 rpm for 1 min
24. Aspirate DNA supernatant and transfer to new tube and freeze.