

プローブベースのqPCR

対応機種/製品	ORA™qPCR Probe Mix, 2X	ORA™qPCR Probe ROX L Mix, 2X	ORA™qPCR Probe ROX H Mix, 2X
	ORA™SEE qPCR Probe Mix, 2X	ORA™SEE qPCR Probe ROX L Mix, 2X	ORA™SEE qPCR Probe ROX H Mix, 2X
Agilent:Mx3000P®,Mx3005P®,Mx4000P®		●	
Analytic Jena:qTOWER	●	●	
BioRad:Opticon®,Opticon®2,Chromo4™,MiniOpticon™,CFX96™,CFX384™	●	●	
BioRad:iCycler®,MyiQ™,iQ™5	●		
BjS:Xpress®	●	●	
Cepheid:SmartCycler®	●	●	
Illumina:Eco	●	●	
Eppendorf:Mastercycler®ep realplex,Mastercycler®realplex 25	●	●	
Fluidigm:BioMark™		●	
Hain Lifescience:FluoroCycler®96	●	●	
Life Technologies:7500,7500 FAST,Viia™7,QuantStudio™12KFlex		●	
Life Technologies:7000,73000,7700,7900,7900HT,7900HT FAST,StepOne™,StepOnePlus™			●
QIAGEN:Rotor-Gene®Q,Rotor-Gene®6000,Rotor-Gene®3000	●	●	
Roche Applied Science:LightCycler®480,LightCycler®96,LightCycler®Nano	●	●	
Takara:Thermal Cycler Dice®	●	●	
Thermo Fisher Scientific:Piko Real®	●	●	
Techne:PrimeQ®,Quantica®	●	●	