

# 色素ベースのRT-qPCR

対応機種/製品	1Step RT qPCR Green ROX L Kit	1Step RT qPCR Green ROX H Kit
Agilent:Mx3000P®,Mx3005P®,Mx4000P®	●	
Analytic Jena:qTOWER	●	
BioRad:Opticon®,Opticon®2,Chromo4™,MiniOpticon™,CFX96™,CFX384™	●	
BioRad:iCycler®,MyiQ™,iQ™5	●	
BjS:Xpress®	●	
Cepheid:SmartCycler®	●	
Illumina:Eco	●	
Eppendorf:Mastercycler®ep realplex,Mastercycler®realplex 25	●	
Fluidigm:BioMark™	●	
Hain Lifescience:FluoroCycler®96	●	
Life Technologies:7500,7500 FAST,Viia™7,QuantStudio™12KFlex	●	
Life Technologies:7000,73000,7700,7900,7900HT,7900HT FAST,StepOne™,StepOnePlus™		●
QIAGEN:Rotor-Gene®Q,Rotor-Gene®6000,Rotor-Gene®3000	●	
Roche Applied Science:LightCycler®480,LightCycler®96,LightCycler®Nano	●	
Takara:Thermal Cycler Dice®	●	
Thermo Fisher Scientific:Piko Real®	●	
Techne:PrimeQ,Quantica®	●	