



Product Focus

ADS™ PCR Cleaning Beads



PCR products, plasmid DNA and Rolling Circle Amplification (RCA) products are the most common templates in Sanger sequencing. Unlike plasmid and RCA templates, which can be directly used for cycle sequencing reactions, PCR products need to be purified before they are used as template for sequencing.

We provide cost-effective ADS™ PCR cleaning beads to purify the PCR products to remove dNTPs and primers in PCR products. Our PCR Cleaning Beads purify PCR product from 100 bp to a few kilobases and produce high-quality and sequencing-ready PCR templates. The product is ideal for high-throughput PCR product clean-up at lower cost compared to spin columns.



ADS™ PCR Cleaning Beads

Catalog Number	Product Name	Unit Size
170001	ADS™ PCR Cleaning Beads	1 ml
170005	ADS™ PCR Cleaning Beads	5 ml