


# Product Focus

## PwrPOP™ Polymers for 3500 Series



PwrPOP™ Polymers (PwrPOP™-4, 6, and 7) for 3500 series are designed specifically for use on AB 3500 or 3500XL. These polymers are packed in disposable pouches and come with new RFID tags for monitoring polymer usage and expiration (see user manual). These polymers are cost-effective alternatives to POP polymers for 3500 series and provide similar sequencing performances. However, 3500 instruments with certain software versions can't use this product; see the user manual for details.



Performance comparison of polymer and running buffer from ADS and Thermo Fisher Scientific on AB 3500XL. Panel 1. AB POP-7 polymer with AB running buffer; panel 2. AB POP-7 polymer with ADS refill running buffer; panel 3. ADS PwrPOP-7 polymer with

### PwrPOP polymers for 3500 series

Catalog numbers	Product name	Sample sizes
340384/340960	PwrPOP™-4 polymer for 3500 Series (No RFID)	384/960
341384/341960	PwrPOP™-4 polymer for 3500 Series (with RFID)	384/960
360384/360960	PwrPOP™-6 polymer for 3500 Series (No RFID)	384/960
361384/361960	PwrPOP™-6 polymer for 3500 Series (with RFID)	384/960
370384/370960	PwrPOP™-7 polymer for 3500 Series (No RFID)	384/960
371384/371960	PwrPOP™-7 polymer for 3500 Series (with RFID)	384/960