

rhIGF-1 LR3

Recombinant Human Insulin-like Growth Factor-1 LR3

Lyophilized powder



Version 6.0 Revision Date: 3.6.2025

Catalog No. HYC007M01

Print Date: 3.6.2025

Introduction

Human Insulin-like growth factor-1 (hIGF-1) is a hormone similar in molecular structure to insulin. It primarily functions through binding to its specific receptor, the insulin-like growth factor 1 receptor (IGF1R). IGF-1 is one of the most potent natural activators of the AKT signaling pathway, which stimulates cell growth and proliferation, and inhibits cell apoptosis. rhIGF-1 LR3 (rhIGF-1 Long Arg3) is an 83 amino acid analog of IGF-1 comprising the complete human IGF-1 sequence with the substitution of an Arg for the Glu at position 3 (hence R3), and a 13 amino acid extension peptide at the N-terminus.

The rhIGF-1 LR3, purified by chromatography from E.coli, then lyophilized with Mannitol, has much higher stability and biological activity. The manufacture capability at hundreds of grams scales can meet your large bioprocessing needs.

Specification

- **Source:** E.coli (*Escherichia coli*)
- **Appearance:** White or off-white loose body
- **Molecular weight:** 8.3 kDa
- **Purity:** ≥90.0% by SDS-PAGE
- **Activity:** ≥ 1.0×10⁶ IU/mg by a cell proliferation assay using Bal B/C 3T3 cell.
- **Endotoxin:** <5.0 EU/mg

Preparation and Storage

- **Reconstitution:** Recommend to reconstitute rhIGF-1 LR3 at 100-200µg/ml concentration in sterile water. Further dilutions can be made in other aqueous buffer.
- **Storage:** Stored at temperature -20°C for at least 36 months.

Please avoid repeated freeze-thaw cycles.

Related Products

| Cat. No. | Product Name |
|-----------|---------------------------------|
| HYC002C01 | OsrHSA Premium, powder |
| HYC002C02 | OsrHSA Premium, liquid |
| HYC002M01 | OsrHSA Advanced, powder |
| HYC002R01 | OsrHSA Precise, powder |
| HYC028M01 | OsrHSA, Fatty acid free, liquid |
| HYC021M01 | OsrhFN, powder |
| HYC044M01 | OsrhTF, powder |
| HYC005M01 | OsrhbFGF, powder |
| HYC026M01 | OsrhVEGF, powder |
| HYC027M01 | OsrhKGF, powder |
| HYC020M01 | rhEGF, powder |
| HYC006M01 | OsrhAAT, powder |
| HYC004M01 | OsrhLF, powder |
| HYC042R01 | OsrhLYZ, powder |
| HYC052R01 | OsrCSA, liquid |
| HYC053R01 | OsrFSA, liquid |

**FOR RESEARCH, LABORATORY AND
MANUFACTURE USE ONLY. NOT INTENDED
FOR DIRECT USE ON HUMANS.**

Wuhan Healthgen Biotechnology Corp.

268 Shendun 5th Road, Wuhan, China 430020

www.oryzogen.net

Email: cjing@oryzogen.com