

**Datasheet 2021/06/22**

**Product Name:** *Anti-Extensin glycoprotein [LM1] Antibody*

**Catalog Number:** *ELD020*

**Size:** *5mL (supernatant)*

**Specifications:**

*Product Type: Antibody*

*Antigen: Extensin glycoprotein*

*Isotype: IgM*

*Clonality: Monoclonal*

*Clone Name: LM1*

*Reactivity: Land plants*

*Immunogen: Polysaccharide*

*Species Immunized: Rat*

*Buffer: Cell Culture Supernatant, 0.05% Sodium Azide*

*Tested Applications: IF (1:10), ELISA (1:10)*

*Storage: +4C (short term), -20C (long term)*

*Shipped: Ambient temperature*

**Provider:** *From the laboratory of [Paul Knox, PhD](#), University of Leeds.*

**References:**

1. *Smallwood, M., Martin, H. & Knox, J.P. An epitope of rice threonine- and hydroxyproline-rich glycoprotein is common to cell wall and hydrophobic plasma-membrane glycoproteins. Planta (1995) 196: 510.*

*If you publish research with this product, please [let us know](#) so we can cite your paper.*

**FOR RESEARCH USE ONLY. NOT INTENDED FOR DIAGNOSTIC or THERAPEUTIC USE.**