

Datasheet 2021/06/22

Product Name: *Anti-Heteromannan [LM22] Antibody*

Catalog Number: *ELD035*

Size: *5mL (supernatant)*

Specifications:

Product Type: Antibody

Antigen: Heteromannan

Isotype: IgM

Clonality: Monoclonal

Clone Name: LM22

Reactivity: Land Plants

Immunogen: Polysaccharide

Species Immunized: Rat

Epitope: Beta-1,4-mannan

Buffer: Cell Culture Supernatant, 0.05% Sodium Azide

Tested Applications: IF 1:10

Storage: -20C

Shipped: Room Temperature

Provider: *From the laboratory of [Paul Knox, PhD](#), University of Leeds.*

References:

1. Marcus SE, Blake AW, Benians TA, et al. Restricted access of proteins to mannan polysaccharides in intact plant cell walls. *Plant J.* 2010;64(2):191-203.

If you publish research with this product, please [let us know](#) so we can cite your paper.

FOR RESEARCH USE ONLY. NOT INTENDED FOR DIAGNOSTIC or THERAPEUTIC USE.