


DNA分子量マーカー選択ガイド


Trans2K[®] Plus DNA Marker (5 µl)



1.5% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM111-xx


Trans2K[®] Plus II DNA Marker (5 µl)



1.5% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM121-xx


Trans5K[®] DNA Marker (5 µl)



1.5% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM141-xx


Trans15K[®] DNA Marker (5 µl)



0.8% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM161-xx


1Kb DNA Ladder (5 µl)



1.0% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM201-xx


1Kb Plus DNA Ladder (5 µl)



1.0% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM211-xx


100bp DNA Ladder (5 µl)



2.0% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM301-xx


Trans2K[®] DNA Marker (5 µl)



1.5% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM101-xx


100bp Plus DNA Ladder (5 µl)



2.0% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM311-xx

100bp Plus II DNA Ladder (5 µl)



2.0% TAE agarose gel stained with GelStain
(loading volume: 5 µl)

Cat#BM321-xx