

# プローブベースのqPCR

対応機種/製品	ORA™qPCR Probe Mix, 2X (品番 : QPP0101、QPP0105)	ORA™qPCR Probe ROX L Mix, 2X (品番 : QPP0201QPP0205)	ORA™qPCR Probe ROX H Mix,2X (品番 : QPP0301、 QPP0305)
	ORA™SEE qPCR Probe Mix, 2X (品番 : QPP0401、QPP0405)	ORA™SEE qPCR Probe ROX L Mix,2X (品番 : QPP0501、QPP0505)	ORA™SEE qPCR Probe ROX H Mix,2X (品番 : QPP0601、 QPP0605)
Agilent:Mx3000P®,Mx3005P®,Mx4000P®		●	
Analytic Jena:qTOWER	●	●	
BioRad:Opticon®,Opticon®2,Chromo4™,MiniOpticon™,CFX96™,CFX384™	●	●	
BioRad:iCycler®,MyiQ™,iQ™5	●		
BjS:Xpress®	●	●	
Cepheid:SmartCycler®	●	●	
Illumina:Eco	●	●	
Eppendorf:Mastercycler®ep realplex,Mastercycler®realpex 25	●	●	
Fluidigm:BioMark™		●	
Hain Lifescience:FluoroCycler®96	●	●	
Life Technologies:7500,7500 FAST,Viia™7,QuantStudio™12KFlex		●	
Life Technologies:7000,73000,7700,7900,7900HT,7900HT FAST,StepOne™,StepOnePlus™			●
QIAGEN:Rotor-Gene®Q,Rotor-Gene®6000,Rotor-Gene®3000	●	●	
Roche Applied Science:LightCycler®480,LightCycler®96,LightCycler®Nano	●	●	
Takara:Thermal Cycler Dice®	●	●	
Thermo Fisher Scientific:Piko Real®	●	●	
Techne:PrimeQ,Quantica®	●	●	