

色素ベースのqPCR

対応機種/製品	ORA™qPCR Green ROX L Mix, 2X	ORA™qPCR Green ROX H Mix, 2X	ORA™qPCR HRM Mix, 2X
	ORA™SEE qPCR Green ROX L Mix, 2X	ORA™SEE qPCR Green ROX H Mix, 2X	
Agilent: Mx3000P®, Mx3005P®, Mx4000P®	●		
Analytic Jena: qTOWER	●		
BioRad: Opticon®, Opticon®2, Chromo4™, MiniOpticon™, CFX96™, CFX384™	●		●
BioRad: iCycler®, MyiQ™, iQ™5	●		
BjS: Xpress®	●		
Cepheid: SmartCycler®	●		
Illumina: Eco	●		●
Eppendorf: Mastercycler®ep realplex, Mastercycler®realplex 25	●		●
Fluidigm: BioMark™	●		
Hain Lifescience: FluoroCycler®96	●		
Life Technologies: 7500, 7500 FAST, Vii™7, QuantStudio™12KFlex	●		●
Life Technologies: 7000, 73000, 7700, 7900, 7900HT, 7900HT FAST, StepOne™, StepOnePlus™		●	●
QIAGEN: Rotor-Gene®Q, Rotor-Gene®6000, Rotor-Gene®3000	●		●
Roche Applied Science: LightCycler®480, LightCycler®96, LightCycler®Nano	●		●
Takara: Thermal Cycler Dice®	●		
Thermo Fisher Scientific: Piko Real®	●		
Techne: PrimeQ, Quantica®	●		