

# 2個買ったら、 3個目無料キャンペーン!!

Elabscience®

※無料品は3点のうち最も安価なものに適用されます。

期間：2019年6月3日～2019年8月30日受注分




## FCM抗体

様々なフローサイトメトリー抗体をラインナップしており、複数の色素から選択可能で、細胞内外の幅広いターゲットを提供しています。




様々なサンプルおよび機器に適合する複数の蛍光標識をラインナップしています。迅速かつ簡便な検出で、ハンドオンタイムは15-20分です。

## アポトーシス検出キット



アポトーシス検出キットはアポトーシス細胞とネクローシス細胞を識別するために開発されました。2種類の蛍光を使用し、早期アポトーシス細胞(アネキシンV+蛍光)、後期アポトーシス細胞またはネクローシス細胞(アネキシンV+蛍光およびDNA+蛍光)、裸核細胞(DNA+蛍光)を検出することが可能です。



左図：Jurkat細胞を1μMカンプトテシンあり(右)またはなし(左)で4時間処理を行なった結果です。



# 一部製品紹介

## FCM抗体

- Anti-Mouse CD19 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD161/NK1.1 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD11c Monoclonal Antibody(各蛍光標識)
- Anti-Human CD83 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD86 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse F4/80 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD16/32 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD29 Monoclonal Antibody(各蛍光標識)
- Anti-Human CD3 Monoclonal Antibody(各蛍光標識)
- Anti-Human CD19 Monoclonal Antibody(各蛍光標識)
- Anti-Human CD16 Monoclonal Antibody(各蛍光標識)
- Anti-Human CD56 Monoclonal Antibody(各蛍光標識)
- Anti-Human CD56 Monoclonal Antibody(各蛍光標識)
- Anti-Human FOXP3 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD62L Monoclonal Antibody(各蛍光標識)
- Anti-Human CD86 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD3 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD132 Monoclonal Antibody(各蛍光標識)
- Anti-Human/Mouse/Rat CD47 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD48 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD54 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD22 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD127/IL-7RA Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD272 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD28 Monoclonal Antibody(各蛍光標識)
- Anti-Mouse CD40 Monoclonal Antibody(各蛍光標識)

他多数ラインナップ

## アポトーシス検出キット

製品	サイズ	品番
Annexin V-FITC/PI Apoptosis Detection Kit	50tests	E-CK-A211-50T
	100tests	E-CK-A211-100T
	200tests	E-CK-A211-200T
Annexin V-FITC/7-ADD Apoptosis Detection Kit	50tests	E-CK-A212-50T
	100tests	E-CK-A212-100T
	200tests	E-CK-A212-200T
Annexin V-AF647/PI Apoptosis Detection Kit	50tests	E-CK-A213-50T
	100tests	E-CK-A213-100T
	200tests	E-CK-A213-200T
Annexin V-AF647/7-AAD Apoptosis Detection Kit	50tests	E-CK-A214-50T
	100tests	E-CK-A214-100T
	200tests	E-CK-A214-200T
Annexin V-PE/7-AAD Apoptosis Detection Kit	50tests	E-CK-A216-50T
	100tests	E-CK-A216-100T
	200tests	E-CK-A216-200T
Annexin V-APC/PI Apoptosis Detection Kit	50tests	E-CK-A217-50T
	100tests	E-CK-A217-100T
	200tests	E-CK-A217-200T
Annexin V-APC/7-AAD Apoptosis Detection Kit	50tests	E-CK-A218-50T
	100tests	E-CK-A218-100T
	200tests	E-CK-A218-200T

フィルジェン 株式会社



【お問い合わせ】 試薬部

TEL : 052-624-4388 FAX : 052-624-4389

メール : biosupport@filgen.jp URL : https://filgen.jp/

代理店

(May.2019)