

キャンペーン期間：2024年3月29日受注分まで

ホルモン・小分子関連 Universal ELISA 割引キャンペーン!!

通常価格

¥148,000-

キャンペーン価格

単品購入時

¥125,800

2点以上購入時

¥96,200



バリデーション済み

スパイクリカバリ、直線性、測定範囲、CV%値、および標準曲線データを検証済みです。

多種の動物由来サンプルに対応

Universal ELISAキットは、ヒト、マウス、非ヒト霊長類、ラット、ブタ、ウシ、その他を含むあらゆる動物種に交差性があります。

高感度

対象のプレコート済みの96ウェル ELISAキットは多くの測定対象項目に対してピコグラムの感度を持っています。

対象製品一覧

さらに、期間中にAssayGenie社製品をご購入の皆様へもれなくメーカーオリジナルのギフト品をプレゼント! ※在庫限り



品番	製品名	サイズ	感度	検出範囲
UNES00055	15d-PGJ2 ELISA	96 assays	1.88ng/mL	3.13-200ng/mL
UNES00011	25-HVD3 (25-Hydroxy Vitamin D3) ELISA	96 assays	1.88ng/mL	3.13-200ng/mL
UNES00025	3-NT (3-Nitrotyrosine) ELISA	96 assays	0.94ng/mL	1.56-100ng/mL
UNES00062	5-HETE (5-Hydroxyeicosatetraenoic Acid) ELISA	96 assays	0.09ng/mL	0.16-10ng/mL
UNES00048	5-HIAA (5-Hydroxyindoleacetic Acid) ELISA	96 assays	0.94ng/mL	1.56-100ng/mL
UNES00035	6-keto-PGF1a (6-keto-prostaglandin F1a) ELISA	96 assays	9.38pg/mL	15.63-1000pg/mL
UNES00065	6-OHDA (6-Hydroxydopamine) ELISA	96 assays	0.09ng/mL	0.16-10ng/mL
UNES00026	8-epi-PGF2 alpha (8-Epi-Prostaglandin F2 Alpha) ELISA	96 assays	9.38pg/mL	15.63-1000pg/mL
UNES00018	8-OHdG (8-Hydroxydeoxyguanosine) ELISA	96 assays	0.94ng/mL	1.56-100ng/mL
UNES00032	Arachidonic Acid ELISA	96 assays	0.94ng/mL	1.56-100ng/mL
UNES00050	Acetylcholine ELISA Kit	96 assays	9.38pg/mL	15.63-1000pg/mL
UNES00027	ADMA (Asymmetrical Dimethylarginine) ELISA	96 assays	9.38ng/mL	15.63-1000ng/mL
UNES00047	Androsterone ELISA	96 assays	0.47ng/mL	0.78-50ng/mL
UNES00045	Aldosterone ELISA	96 assays	18.75pg/mL	31.25-2000pg/mL

品番	製品名	サイズ	感度	検出範囲
UNES00049	Bilirubin ELISA Kit	96 assays	0.19Mug/mL	0.31-20Mug/mL
UNES00060	Dihydrobiopterin (BH2) ELISA	96 assays	18.75pg/mL	31.25-2000pg/mL
UNES00036	cAMP (Cyclic AMP) ELISA	96 assays	0.94ng/mL	1.56-100ng/mL
UNES00044	Carnosine ELISA	96 assays	9.38ng/mL	15.63-1000ng/mL
UNES00052	cGMP (Cyclic GMP) ELISA	96 assays	0.47pmol/mL	0.78-50pmol/mL
UNES00038	Creatinine ELISA	96 assays	0.75Mug/mL	1.25-80Mug/mL
UNES00029	Dopamine ELISA	96 assays	18.75pg/mL	31.25-2000pg/mL
UNES00021	DHT (Dihydrotestosterone) ELISA	96 assays	23.44pg/mL	39.06-2500pg/mL
UNES00012	DHVD3 (1,25-Dihydroxyvitamin D3) ELISA	96 assays	4.69pg/mL	7.81-500pg/mL
UNES00028	Epinephrine (Adrenaline) ELISA	96 assays	18.75pg/mL	31.25-2000pg/mL
UNES00006	Folic Acid (Vitamin B9) ELISA	96 assays	0.94ng/mL	1.56-100ng/mL
UNES00046	GnRH (Gonadotropin-Releasing Hormone) ELISA	96 assays	9.38pg/mL	15.63-1000pg/mL
UNES00017	GSH (Glutathione) ELISA	96 assays	0.94Mug/mL	1.56-100Mug/mL
UNES00022	Histamine ELISA	96 assays	0.94ng/mL	1.56-100ng/mL
UNES00039	Inositol Triphosphate (IP3) ELISA	96 assays	9.38pg/mL	15.63-1000pg/mL
UNES00041	Leukotriene B4 (LTB4) ELISA	96 assays	9.38pg/mL	15.63-1000pg/mL
UNES00061	Leukotriene E4 (LTE4) ELISA	96 assays	18.75pg/mL	31.25-2000pg/mL
UNES00034	Lipoxin A4 (LXA4) ELISA	96 assays	0.47ng/mL	0.78-50ng/mL
UNES00040	MDA (Malondialdehyde) ELISA	96 assays	18.75ng/mL	31.25-2000ng/mL
UNES00015	Methionine-Enkephalin (MEK) ELISA	96 assays	75pg/mL	125-8000pg/mL
UNES00030	Noradrenaline (Norepinephrine) ELISA	96 assays	0.19ng/mL	0.31-20ng/mL
UNES00019	Oxytocin ELISA	96 assays	9.38pg/mL	15.63-1000pg/mL
UNES00066	Progesterone ELISA Kit	96 assays	0.15ng/mL	0.31-20ng/mL
UNES00043	Prostaglandin D2 (PGD2) ELISA	96 assays	9.38pg/mL	15.63-1000pg/mL
UNES00033	Prostaglandin E1 (PGE1) ELISA	96 assays	46.88pg/mL	78.13-5000pg/mL
UNES00024	Prostaglandin E2 (PGE2) ELISA	96 assays	18.75pg/mL	31.25-2000pg/mL
UNES00016	Prostacyclin ELISA	96 assays	18.75pg/mL	31.25-2000pg/mL
UNES00054	Pregnenolone ELISA	96 assays	0.47ng/mL	0.78-50ng/mL
UNES00057	Pentosidine (PTD) ELISA	96 assays	0.47ng/mL	0.78-50ng/mL
UNES00023	Serotonin (5-HT) ELISA	96 assays	9.38ng/mL	15.63-1000ng/mL
UNES00063	STAT-3 ELISA Kit	96 assays	0.19ng/mL	0.31-20ng/mL
UNES00037	Thromboxane A2 (TXA2) ELISA	96 assays	9.38pg/mL	15.63-1000pg/mL
UNES00005	Vitamin B6 ELISA	96 assays	0.94ng/mL	1.56-100ng/mL
UNES00009	Vitamin D ELISA	96 assays	3.75ng/mL	6.25-400ng/mL
UNES00010	Vitamin D3 ELISA	96 assays	0.94ng/mL	1.56-100ng/mL
UNES00014	Vitamin K1 ELISA	96 assays	9.38pg/mL	15.63-1000pg/mL

フィルジェン 株式会社



【お問い合わせ】 試薬機器部

TEL : 052-624-4388 FAX : 052-624-4389

メール : biosupport@filgen.jp URL : https://filgen.jp/

代理店

(DEC.2023)